



MGUARD

MGUARD FACESHIELD

MGuard Faceshield

Item	Face shield
Color	Transparent
Material	PET
Size	32 x 22cm
Style	Full Length
Layer	3-ply(PET shield +2-ply PE film)
Component	Plastic Shield, foam, elastic straps
Minimum Order Quantity	250 PCS
Package	1 per LDPE bag, 250 bags per carton



Description - Faceshield With Foam Headband & Elastic Straps

- Foam headband forms a seal between forehead and shield
- Shield has an antifog coating
- Crystal clear viewing
- Good peripheral coverage
- Shield sits off the face allowing the use of glasses
- Individually bagged
- Listed with the TGA

